


Product Datasheet

Estrogen Receptor beta antibody (orb1473715)

Description	Mouse monoclonal antibody to Estrogen Receptor beta
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Estrogen Receptor beta. The exact sequence is proprietary.
Target	ESR2
Preservatives	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q92731
Entrez	2100
Dilution Range	WB (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.



Western blot analysis of Estrogen Recept...