

Product Datasheet

Acetylated Lysine antibody (orb147372)

Description	Mouse monoclonal to Acetylated Lysine (ATTO-488). Post-translational modifications of proteins play ...
Species/Host	Mouse
Reactivity	All
Conjugation	ATTO-488
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	Acetylated KLH
Target	Acetylated Lysine
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg of SMC-153 was sufficient to detect acetylated chicken erythrocyte histones (sodium butyrate-treated) using 20 µg total protein, on western blot by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB (1:1000), IHC (1:100)
Expiration Date	12 months from date of receipt.