



Product Datasheet

CD366 antibody (orb1473721)

Description	Mouse monoclonal antibody to CD366
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant fusion protein of human CD366. The exact sequence is proprietary.
Target	HAVCR2
Preservatives	Mouse IgG2b kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG2b kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q8TDQ0
Entrez	84868
Dilution Range	WB (1/500 - 1/2000), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of CD366 expression...



Flow cytometric analysis of U2OS cells u...