

Product Datasheet

Acetylated Lysine antibody [Out of stock] (orb147374)

Description

Mouse monoclonal to Acetylated

modifications of proteins play ...

Species/Host

Mouse

Reactivity

All

Conjugation

ATTO-594

Tested Applications

ELISA, ICC, IHC, IP, WB

Immunogen

Acetylated KLH

Target

Acetylated Lysine

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg of SMC-153 was sufficient to detect acetylated chicken erythrocyte histones (sodium butyrate-treated) using 20 µg total protein, on western blot by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

Dilution Range

WB (1:1000), IHC (1:100)

Expiration Date

12 months from date of receipt.