


Product Datasheet

PAX8 antibody (orb1473802)

Description	Mouse monoclonal antibody to PAX8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PAX8. The exact sequence is proprietary.
Target	PAX8
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q06710
Entrez	7849
Dilution Range	WB (1/1000 - 1/2000)
Expiration Date	12 months from date of receipt.



Western blot analysis of PAX8 expression...