



Product Datasheet

Alpha-actin-1 antibody (orb1473825)


Description	Mouse monoclonal antibody to Alpha-actin-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	Recombinant fusion protein of human Alpha-actin-1. The exact sequence is proprietary.
Target	ACTA1
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P68133
Entrez	58
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of Alpha-actin-1 e...



Immunohistochemical analysis of Alpha-ac...



Flow cytometric analysis of A549 cells u...