



Product Datasheet

ALDH3A2 antibody (orb1473868)

Description	Mouse monoclonal antibody to ALDH3A2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant fusion protein of human ALDH3A2. The exact sequence is proprietary.
Target	ALDH3A2
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P51648
Entrez	224
Dilution Range	WB (1/1000 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of ALDH3A2 express...



Immunohistochemical analysis of ALDH3A2

...