



Product Datasheet

FHIT antibody (orb1473881)

Description	Mouse monoclonal antibody to FHIT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant fusion protein of human FHIT. The exact sequence is proprietary.
Target	FHIT
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P49789
Entrez	2272
Dilution Range	WB (1/500 - 1/2000), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of FHIT expression...



Flow cytometric analysis of HepG2 cells ...