



Product Datasheet

CD49f antibody (orb1473946)


Description	Mouse monoclonal antibody to CD49f
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IH, WB
Immunogen	Recombinant fusion protein of human CD49f. The exact sequence is proprietary.
Target	ITGA6
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P23229
Entrez	3655
Dilution Range	WB (1/1000 - 1/5000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of CD49f expressio...



Immunohistochemical analysis of CD49f st...



Immunofluorescent analysis of CD49f stai...