



Product Datasheet

CD107a antibody (orb1474018)


Description	Mouse monoclonal antibody to CD107a
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human CD107a. The exact sequence is proprietary.
Target	LAMP1
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P11279
Entrez	3916
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of CD107a expressi...



Immunohistochemical analysis of CD107a S...



Immunofluorescent analysis of CD107a sta...