



Product Datasheet

MAP2 antibody (orb1474020)


Description	Mouse monoclonal antibody to MAP2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human MAP2. The exact sequence is proprietary.
Target	MAP2
Preservatives	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P11137
Entrez	4133
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of MAP2 expression...



Immunohistochemical analysis of MAP2 staining...



Immunofluorescent analysis of MAP2 staining...