



Product Datasheet

CD10 antibody (orb1474040)


Description	Mouse monoclonal antibody to CD10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human CD10. The exact sequence is proprietary.
Target	MME
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P08473
Entrez	4311
Dilution Range	WB (1/500 - 1/2000), IH (1/100 - 1/400), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of CD10 expression...



Immunohistochemical analysis of CD10 sta...



Flow cytometric analysis of Jurkat cells...