



Product Datasheet

HSP72 antibody (orb1474047)

Description	Mouse monoclonal antibody to HSP72
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human HSP72. The exact sequence is proprietary.
Target	HSPA1A; HSPA1B
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P0DMV8 , P0DMV9
Entrez	3304 , 3303
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of HSP72 expression...



Immunofluorescent analysis of HSP72 staining...