



Product Datasheet

RET antibody (orb1474050)


Description	Mouse monoclonal antibody to RET
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human RET. The exact sequence is proprietary.
Target	RET
Preservatives	Mouse IgM kappa. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgM kappa. Supplied in crude ascites with 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P07949
Entrez	5979
Dilution Range	WB (1/500 - 1/5000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of RET expression ...



Immunohistochemical analysis of RET stai...



Immunofluorescent analysis of RET staini...