

Product Datasheet

CD29 antibody (orb1474071)

Description	Mouse monoclonal antibody to CD29
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CD29. The exact sequence is proprietary.
Target	ITGB1
Preservatives	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P05556
Entrez	3688
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.

