


Product Datasheet

IGHD antibody (orb1474109)

Description	Mouse monoclonal antibody to IGH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human IGH. The exact sequence is proprietary.
Target	IGH
Preservatives	Mouse IgM. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgM. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P01880
Dilution Range	WB (1/100 - 1/500)
Expiration Date	12 months from date of receipt.



Western blot analysis of IGH expression...