



Product Datasheet

PGK1 antibody (orb1474128)


Description	Mouse monoclonal antibody to PGK1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PGK1. The exact sequence is proprietary.
Target	PGK1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P09411, P00558
Entrez	18655, 5230
Dilution Range	WB (1/500 - 1/2000), IH (1/50 - 1/100), FC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of PGK1 expression...



Immunohistochemical analysis of PGK1 staining...



Overlay histogram showing Jurkat cells S...