



Product Datasheet

EEF1E1 antibody (orb1474159)


Description	Mouse monoclonal antibody to EEF1E1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant fusion protein of human EEF1E1. The exact sequence is proprietary.
Target	EEF1E1
Preservatives	Mouse IgG2b kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Mouse IgG2b kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	O43324 , Q9D1M4
Entrez	66143 , 9521
Dilution Range	WB (1/500 - 1/2000), IF/IC (1/10 - 1/50), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of EEF1E1 expressi...



Immunofluorescent analysis of EEF1E1 sta...



Flow cytometric analysis of HeLa cells u...