



Product Datasheet

SIRT6 antibody (orb1474222)

Description	Mouse monoclonal antibody to SIRT6
Species/Host	Mouse
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Purified recombinant human SIRT6 protein expressed in E.coli.
Target	SIRT6
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q8N6T7 , P59941
Entrez	50721 , 51548
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/100), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of SIRT6 expression...



Immunofluorescent analysis of SIRT6 staining...