



Product Datasheet

Metadherin antibody (orb1474245)


Description	Mouse monoclonal antibody to Metadherin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	Purified recombinant fragment of human Metadherin expressed in E. Coli.
Target	MTDH
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q86UE4
Entrez	92140
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), FC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of Metadherin expr...



Immunohistochemical analysis of Metadher...



Flow cytometric analysis of Metadherin i...