



Product Datasheet

WDR77 antibody (orb1474311)

Description	Mouse monoclonal antibody to WDR77
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Purified recombinant human WDR77 protein fragments expressed in E.coli.
Target	WDR77
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q99J09 , Q9BQA1 , Q4QR85
Entrez	310769 , 70465 , 79084
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of WDR77 expression...



Immunofluorescent analysis of WDR77 staining...