



Product Datasheet

PDC-E2 antibody (orb1474342)

Description	Mouse monoclonal antibody to PDC-E2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Purified recombinant human Pyruvate Dehydrogenase E2 protein fragments expressed in E.coli.
Target	DLAT
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P10515 , Q8BMF4
Entrez	1737 , 235339
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/100), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of PDC-E2 expressi...



Immunofluorescent analysis of PDC-E2 sta...