


## Product Datasheet

### CKMT1A antibody (orb1474350)

|                            |                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CKMT1A                                                              |
| <b>Species/Host</b>        | Mouse                                                                                            |
| <b>Reactivity</b>          | Human, Mouse                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                     |
| <b>Tested Applications</b> | WB                                                                                               |
| <b>Immunogen</b>           | Purified recombinant human CKMT1 protein fragments expressed in E.coli.                          |
| <b>Target</b>              | CKMT1A                                                                                           |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.             |
| <b>Form/Appearance</b>     | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.             |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                            |
| <b>Clonality</b>           | Monoclonal                                                                                       |
| <b>Source</b>              | Mouse                                                                                            |
| <b>Uniprot ID</b>          | <b>P12532, P30275</b>                                                                            |
| <b>Entrez</b>              | <b>1159, 12716</b>                                                                               |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000)                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                  |



Western blot analysis of CKMT1A expressi...