



Product Datasheet

TBL1XR1 antibody (orb1474382)


Description	Mouse monoclonal antibody to TBL1XR1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Purified recombinant human TBLR1 protein fragments expressed in E.coli.
Target	TBL1XR1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q8BHJ5 , Q9BZK7
Entrez	79718 , 81004
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of TBL1XR1 express...



Immunohistochemical analysis of TBL1XR1
...



Immunofluorescent analysis of TBL1XR1
st...