


Product Datasheet

LC3-I/II antibody (orb1474437)

Description	Rabbit monoclonal antibody to LC3-I/II
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human MAP1LC3A
Target	MAP1LC3A
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	Q6XVN8 , Q9H492 , Q91VR7
Entrez	66734 , 362245 , 84557
Dilution Range	WB (1/500 - 1/1000), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of LC3-I/II expres...