

Product Datasheet

Filamin A/B (phospho-S2152/2107) antibody (orb1474464)

Description

Rabbit monoclonal antibody to Filamin A/B

kDa A B C D

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IF, IH, WB

Applications

Western blot analysis
of Filamin A/B (Ph...

Immunogen


KLH-conjugated synthetic phosphopeptide
corresponding to residues surrounding
S2152/2107 of human Filamin A/B protein.
The exact sequence is proprietary.

Target

FLNB

Preservatives


Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M
NaCl, 50% Glycerol, 0.01% Sodium azide
and 0.05% rAlbumin.



Immunohistochemical
analysis of Filamin ...

Form/Appearance

Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M
NaCl, 50% Glycerol, 0.01% Sodium azide
and 0.05% rAlbumin.



Immunofluorescent
analysis of Filamin
A/...

Storage

Shipped at 4°C. Upon delivery aliquot and
store at -20°C for one year. Avoid
freeze/thaw cycles.

Note

For research use only

Clonality

Monoclonal

Source

Rabbit

Uniprot ID

075369

Entrez

2317

Dilution Range

WB (1/500 - 1/1000), IH (1/50 - 1/100), IF/IC
(1/50 - 1/100)

Expiration Date

12 months from date of receipt.