



Product Datasheet

Histone H3 (AcK27) antibody (orb1474466)


Description	Rabbit monoclonal antibody to Histone H3 (Acetyl-K27)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K27 of human Histone H3 protein. The exact sequence is proprietary.
Target	H3C1
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	P68431
Entrez	8350
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), IF/IC (1/50 - 1/100), IP (1/10 - 1/50), ChIP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (Ace...



Immunohistochemical analysis of Histone ...



Immunofluorescent analysis of Histone H3...