



Product Datasheet

HSP90 alpha/beta antibody (orb1474524)


Description	Rabbit monoclonal antibody to HSP90 alpha/beta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	Recombinant protein of human Hsp90 beta
Target	HSP90AA1; HSP90AB1
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	P07901 , P82995 , P34058 , P11499 , P07900 , P08238
Entrez	3326 , 15516 , 3320 , 15519 , 301252 , 299331
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), IF/IC (1/50 - 1/100), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of HSP90 alpha/bet...



Immunohistochemical analysis of HSP90 al...



Immunofluorescent analysis of HSP90 alph...