



Product Datasheet

GSK3 beta (phospho-S9) antibody (orb1474580)


Description	Rabbit monoclonal antibody to GSK3 beta (phospho-S9)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S9 of human GSK3 beta protein. The exact sequence is proprietary.
Target	GSK3B
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	P49841, P18266, Q9WV60
Entrez	84027, 2932, 56637
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of GSK3 beta (Phos...



Immunohistochemical analysis of GSK3 bet...



Immunofluorescent analysis of GSK3 beta ...