



Product Datasheet

GSK3 beta antibody (orb1474591)

Description	Rabbit monoclonal antibody to GSK3 beta
Species/Host	Rabbit
Reactivity	Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	A synthetic peptide of human GSK3 beta
Target	GSK3B
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Form/Appearance	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Rabbit
Uniprot ID	P49841 , P18266
Entrez	84027 , 2932
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100)
Expiration Date	12 months from date of receipt.



Western blot analysis of GSK3 beta expression.



Immunohistochemical analysis of GSK3 beta.