



Product Datasheet

TRAF6 antibody (orb1474642)

Description	Rabbit polyclonal antibody to TRAF6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TRAF6. The exact sequence is proprietary.
Target	TRAF6
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P70196, B5DF45, Q9Y4K3
Entrez	7189, 22034, 311245
Dilution Range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of TRAF6 expression...



Immunofluorescent analysis of TRAF6 staining...