



Product Datasheet

MUC2 antibody (orb1474645)


Description	Rabbit polyclonal antibody to MUC2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human MUC2. The exact sequence is proprietary.
Target	MUC2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q80Z19 , Q62635 , Q02817
Entrez	17831 , 24572 , 4583
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), IF/IC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of MUC2 expression...



Immunohistochemical analysis of MUC2 sta...



Immunofluorescent analysis of MUC2 stain...