


Product Datasheet

Histone H4 (AcK91) antibody (orb1474655)

Description	Rabbit polyclonal antibody to Histone H4 (Acetyl-K91)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K91 of human Histone H4 protein. The exact sequence is proprietary.
Target	HIST1H4A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P62805
Entrez	8361, 121504, 8367, 8365, 8359, 8364, 8366, 8360, 8362, 8370, 8368, 8294, 8363, 554313
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H4 (Ace...



Direct ELISA antibody dose-response curv...