






Product Datasheet

TST Antibody (orb1474807)

Description	TST Antibody	
Species/Host	Rabbit	  IHC analysis of TST using anti-TST antib...
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, ICC, IF, IHC, WB	
Immunogen	E.coli-derived human TST recombinant protein (Position: M1-D220).	 IHC analysis of TST using anti-TST antib...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	 IHC analysis of TST using anti-TST antib...
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	34 kDa	 IHC analysis of TST using anti-TST antib...
Uniprot ID	Q16762	
Dilution Range	Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Direct ELISA, 0.1-0.5 µg/ml, Human	
Expiration Date	12 months from date of receipt.	