




Product Datasheet

SNRPA1 Antibody (orb1474822)

Description	SNRPA1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, ICC, IF, IHC, WB	IHC analysis of SNRPA1 using anti-SNRPA1...
Immunogen	E.coli-derived human SNRPA1 recombinant protein (Position: M1-E233).	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	IHC analysis of SNRPA1 using anti-SNRPA1...
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Isotype	Rabbit IgG	IHC analysis of SNRPA1 using anti-SNRPA1...
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	P09661	
Dilution Range	Western blot, 0.1-0.25 µg/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human, Mouse, Rat Flow Cytometry, 1-3 µg/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human	
Expiration Date	12 months from date of receipt.	