

Product Datasheet

SNRNP27 Antibody (orb1474826)

Description

SNRNP27 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC, WB

Applications

Immunogen

E.coli-derived human SNRNP27 recombinant protein (Position: E100-A155).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

23 kDa

Uniprot ID



[Q8WVK2](https://www.uniprot.org/entry/Q8WVK2)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human
Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of SNRNP27 using anti-SNRNP...



IHC analysis of SNRNP27 using anti-SNRNP...



IHC analysis of SNRNP27 using anti-SNRNP...