

Product Datasheet

Desmoglein 2/DSG2 Antibody (monoclonal, 2B4D1) (orb1474876)

Description

Desmoglein 2/DSG2 Antibody (monoclonal, 2B4D1)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E.coli-derived human Desmoglein 2/DSG2 recombinant protein (Position: L24-E1020).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

IgG1

Clonality

Monoclonal

MW


160 kDa

Uniprot ID
[Q14126](#)
Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section),
2-5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis
of
Desmoglein
2/DSG2
using ...



IHC analysis
of
Desmoglein
2/DSG2
using ...



IHC analysis
of
Desmoglein
2/DSG2
using ...