

Product Datasheet

Golgin 97/GOLGA1 Antibody (monoclonal, 8E4H1) (orb1474879)

Description

Golgin 97/GOLGA1 Antibody (monoclonal, 8E4H1)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IF, WB

Applications

Immunogen

E.coli-derived human Golgin 97/GOLGA1 recombinant protein (Position: M1-K752).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

IgG1

Clonality

Monoclonal

MW

97 kDa

Uniprot ID


[Q92805](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



IF analysis of Golgin 97/GOLGA1 using an...



Western blot analysis of Golgin 97/GOLGA...



Flow Cytometry analysis of A549 cells us...