


Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, FITC Conjugated (orb1474924)

Description	Goat Anti-Rabbit IgG (H+L) Secondary Antibody, FITC Conjugated
Species/Host	Goat
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Rabbit IgG (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



IF
analysis
of rat
brain
tissues
using r...