

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Secondary Antibody, Biotin Conjugated (orb1474935)

Description

Goat Anti-Rabbit IgG (H+L) Secondary



Species/Host

Goat

Reactivity

Rabbit

Conjugation

Biotin

Tested Applications

ELISA, IHC

Immunogen

Rabbit IgG (whole molecule).

Concentration

1mg/ml

Storage

At -20°C for one year from date of receipt.
Avoid repeated freezing and thawing.

Note


For research use only

Clonality


Polyclonal

Expiration Date


12 months from date of receipt.



IHC
analysis
of mouse
testis
tissues
usi...



IHC
analysis
of human
intestinal
cancer ...



IHC
analysis
of mouse
liver
tissues
usin...