



Product Datasheet

Goat Anti-Mouse IgG (H+L) Secondary Antibody, Biotin Conjugated (orb1474936)

Description	Goat Anti-Mouse IgG (H+L) Secondary Antibody, Biotin Conjugated
Species/Host	Goat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Mouse IgG (whole molecule).
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



IHC analysis of rat brain tissues using ...



IHC analysis of human intestinal cancer ...