

Product Datasheet

Human SPINT1 ELISA Kit (orb1475021)

Description	Human Kunitz-type Protease Inhibitor 1 (SPINT1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O43278
Sensitivity	0.061ng/mL
Entrez	6692
Expiration Date	Please enquire.