

Product Datasheet

Mouse LAB7-2 ELISA Kit (orb1475035)

Description	Mouse B-Lymphocyte Activation Antigen B7-2 (LAB7-2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42082
Sensitivity	0.058ng/mL
Entrez	12524
Expiration Date	Please enquire.