

Product Datasheet

HSP27 antibody (orb147518)

Description

Mouse monoclonal to Hsp27 (FITC). Hsp27s

	Hsp27s belong...
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	ELISA, FC, ICC, IHC, IP, WB
Immunogen	Full length human HSP27
Target	HSP27
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.5 µg/ml of SMC-161 was sufficient for detection of HSP27 in 10 µg of HeLa lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2b Kappa
Clonality	Monoclonal
MW	27kDa
Uniprot ID	P04792
NCBI	NP_001531.1
Entrez	3315
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.