

Product Datasheet

Human IGFBP-5 Protein (Biotin) (orb1475199)

Description	Human IGFBP-5 (C-6His-Avi) conjugated to Biotin
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Conjugation	Biotin
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 300mM NaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 300mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His-Avi
Note	For research use only
Protein Sequence	Leu21-Glu272
Source	Human Cells
Uniprot ID	P24593
Expiration Date	6 months from date of receipt.