



Product Datasheet

NFκB-p65 (Phospho-Thr254) Antibody (orb14752)


Description	NFκB-p65 (Phospho-Thr254) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q04206
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal Immunofluorescent analysis of H...



Western blot analysis of NFκB-p65(Phospho-Thr254)