

Product Datasheet

Human NTNG1 Protein (Biotin) (orb1475201)

Description	Human NTNG1 (C-Avi-6His) conjugated to Biotin
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Conjugation	Biotin
Preservatives	Lyophilized from a 0.2 μm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Tag	C-Avi-6His
Note	For research use only
Protein Sequence	His29-Ser409
Source	Human Cells
Uniprot ID	Q9Y2I2
Expiration Date	6 months from date of receipt.