

Product Datasheet

Human PVRIG Protein (Biotin) (orb1475202)

Description	Human PVRIG (C-Fc-Avi) conjugated to Biotin
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Conjugation	Biotin
Preservatives	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Tag	C-Fc-Avi
Note	For research use only
Protein Sequence	Thr41-Asp171
Source	Human Cells
Uniprot ID	Q6DKI7
Expiration Date	6 months from date of receipt.