

Product Datasheet

Human LAIR2 Protein (Biotin) (orb1475233)

Description	Human LAIR2 (C-Avi-6His) conjugated to Biotin
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Conjugation	Biotin
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Avi-6His
Note	For research use only
Protein Sequence	Gln22-Pro152
Source	Human Cells
Uniprot ID	Q6ISS4
Expiration Date	6 months from date of receipt.