

Product Datasheet

2019-nCoV Papain-Like Protease Protein (orb1475260)

Description	Recombinant 2019-nCoV Papain-Like Protease
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-Hcl, 10 mM 2-Mercaptoethanol, 20%glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-Hcl, 10 mM 2-Mercaptoethanol, 20%glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Glu1564-Lys1878
Source	E.coli
Uniprot ID	QHD43415.1
Expiration Date	6 months from date of receipt.